

## Campylobacter sp. Selective Agar (CB)

In vitro diagnostic medical aid.

Only for professional use.

Sterile

**Purpose:** The product is intended for the selective cultivation and isolation of Campylobacter sp. at 42 °C and Arcobacter sp. At 25 °C.

**Package contents:** 18 ml in a Petri dish

**Composition:** Peptone, heart infusion, gelatin, starch, calcium hydroxide, buffer, agar, cefoperazone, vancomycin, sheep blood, demineralized water

**Storage:** Store in a cool place at +4 to +6 °C.

**Usability:** Use within 24 hours after opening the primary packaging.

**Procedure:** Inoculate the clinical material aseptically on the medium surface.

Incubation: under optimal conditions – 24 to 48 hours in an oxygen-depleted atmosphere enriched with carbon dioxide at temperatures of: C. fetus ssp. fetus at 25 °C and 37 °C, C. jejuni/coli at 37 °C and 42 °C, C. fetus spp. venerealis at 25 °C and 37 °C.

**Information for users:** Do not use the medium if it shows signs of microbial contamination, color change, drying, cracking or other signs of damage. All specimens and inoculated products must be considered as potentially infectious and safety measures must be followed in accordance with applicable regulations when handling them.

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

**Instructions in case of a damaged packaging:** In case of damage to the packaging, do not use and dispose of the product.

### Quality control:

Type of inspection	Inoculated culture	Result
positive check	Campylobacter jejuni CCM7212	culture growth
negative check	Escherichia coli ATCC 25922 / CCM 3954	growth inhibition

### Manufacturer

EnviroLab s.r.o.  
Galvaniho 7  
821 04 Bratislava  
Slovak Republic

### Place of business:

Obchodná ulica 7  
946 51 Nesvady  
Slovak Republic