
Christensen's Agar – slant urea (CHM)

In vitro diagnostic medical aid.

Only for professional use.

Sterile

Purpose: The product is used to detect production of urease.

Package contents: 4 ml in a test tube

Composition: Agar base, urea, demineralized water

Storage: Store in a cool place at +4 to +6 °C.

Usability: Use within 24 hours after opening the primary packaging.

Procedure: Inoculate the clinical material aseptically on the medium surface.

Incubation: under optimal conditions – 24 hours at 35 °C.

Information for users: Do not use the medium if it shows signs of microbial contamination, color change, drying, cracking or other signs of damage. All specimens and inoculated products must be considered as potentially infectious and safety measures must be followed in accordance with applicable regulations when handling them.

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

| Type of inspection | Inoculated culture | Result |
|--------------------|----------------------------|----------------|
| positive check | Trichophyton interdigitale | cyclamen color |
| negative check | Trichophyton rubrum | yellow color |

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic