

Triple Sugar Iron Agar (HAJN)

In vitro diagnostic medical aid.

Only for professional use.

Sterile

Purpose: The product is designed for the detection of basic biochemical properties of gram-negative bacteria.

Package contents: 5 or 7 ml in a test tube

Composition: Meat extract, peptone, yeast extract, lactose, glucose, sucrose, sodium chloride, phenol red, sodium thiosulfate, ferricammonium citrate, agar, demineralized water

Storage: Store in a cool place at +4 to +6 °C.

Usability: Use within 24 hours after opening the primary packaging.

Procedure: Inoculate the clinical material aseptically on the medium surface.

Incubation: under optimal conditions – 24 hours at 35 °C.

Information for users: Do not use the medium if it shows signs of microbial contamination, color change, drying, cracking or other signs of damage. All specimens and inoculated products must be considered as potentially infectious and safety measures must be followed in accordance with applicable regulations when handling them.

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

Type of inspection	Inoculated culture	Result
positive check	Escherichia coli CCM3954 Proteus mirabilis CCM1944 Shigella flexneri CCM4422 Pseudomonas aeruginosa CCM3955	growth with gas formation, fermented glycolides
negative check	uninoculated medium	without a change

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic