

---

**Selective Neisseria sp. Agar with Vancomycin (NSA+ATB)**


---

In vitro diagnostic medical aid.

**Only for professional use.**

**Sterile**

---

**Purpose:** The product is designed for the selective isolation of *Neisseria gonorrhoeae* and *Neisseria meningitidis*.

**Package contents:** 18 ml in a Petri dish

**Composition:** Peptone, heart infusion, gelatin, starch, calcium hydroxide, buffer, agar, colistin sulfate, vancomycin hydrochloride, nystatin, trimethoprim, sheep blood, demineralized water

**Storage:** Store in a cool place at +4 to +6 °C.

**Usability:** Use within 24 hours after opening the primary packaging.

**Procedure:** Inoculate the clinical material aseptically on the medium surface.

Incubation: under optimal conditions – 24 hours at 35 °C.

**Information for users:** Do not use the medium if it shows signs of microbial contamination, color change, drying, cracking or other signs of damage. All specimens and inoculated products must be considered as potentially infectious and safety measures must be followed in accordance with applicable regulations when handling them.

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

**Instructions in case of a damaged packaging:** In case of damage to the packaging, do not use and dispose of the product.

**Quality control:**

Type of inspection	Inoculated culture	Result
positive check	<i>Neisseria gonorrhoeae</i> ATCC19424/CCM4500	culture growth
negative check	<i>Enterococcus faecalis</i> CCM4224 <i>Candida albicans</i> CCM8215	growth inhibition

**Manufacturer**

EnviroLab s.r.o.  
Galvaniho 7  
821 04 Bratislava  
Slovak Republic

**Place of business:**

Obchodná ulica 7  
946 51 Nesvady  
Slovak Republic