

Sabouraud Agar with 4% Glucose (SAB)

In vitro diagnostic medical aid.

Only for professional use.

Sterile

Purpose: The product is intended for the isolation and cultivation of microscopic fungi from a clinical material.

Package contents: 18 ml in a Petri dish

Composition: peptone, glucose, agar, demineralized water

Storage: Store in a cool place at +4 to +6 °C.

Usability: Use within 24 hours after opening the primary packaging.

Procedure: Inoculate the sample directly onto the surface of the medium; incubate at 25 °C for 5-7 days. After cultivations, identify the colonies using a microscopic specimen, biochemical tests, and based on the assimilation and fermentation properties.

Information for users: Do not use the medium if it shows signs of microbial contamination, color change, drying, cracking or other signs of damage. All specimens and inoculated products must be considered as potentially infectious and safety measures must be followed in accordance with applicable regulations when handling them.

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

Type of inspection	Inoculated culture	Result
positive check	Candida albicans CCM8215	growth, creamy white smooth convex colonies
	Aspergillus niger CCM8189	growth, larger colonies with black sporulation
negative check	Staphylococcus aureus CCM3953	no growth
	Escherichia coli CCM3954	

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic