
Müller-Hinton Agar (MHB)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed to prepare a bacterial suspension to detect the minimum inhibitory concentration of antibiotics.

Package contents: 8 ml and 10 ml in a test tube.

Composition: Dehydrated meat infusion, casein hydrolyzate, starch, demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 24 hours after opening the primary packaging.

Procedure: Inoculate the clinical material aseptically in the growth medium..

Incubation: under optimal conditions – 24 hours at 35 °C.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

Type of	Inoculated culture	Result
Positive inspection	<i>Staphylococcus aureus</i> ATCC25923 / CCM3953 <i>Escherichia coli</i> ATCC 25922 /CCM 3954 <i>Pseudomonas aeruginosa</i> CCM3955 <i>Enterococcus faecalis</i> CCM4224	culture growth
Negative check	uninoculated medium	without a change

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic