
Glucose 20% (GLU20)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed for the transport of samples.

Package contents: 1 ml in a test tube

Composition: Glucose, demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 24 hours after opening the primary packaging.

Procedure: Inoculate the clinical material aseptically into the solution.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic