

0.1% Carbol-Fuchsin Solution (KARB)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed for Gram staining of microscopic preparations.

Package contents: 1000 ml solution in a bottle

Composition: ethanol, fuchsin (C.I. 42510), phenol

demineralized water

Storage: Store in a dark bottle below +25 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users:

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Warning

H226 Flammable liquid and vapor
H302 Harmful if swallowed
H312 Harmful in contact with skin
H332 Harmful if inhaled

P210	Keep away from heat / sparks / open flame / hot surfaces. Do not smoke.
P233	Keep the container tightly closed.
P240	Ground/fix the container and the filler.
P241	Use an electrical / ventilation / lighting device designed for an explosive atmosphere.
P242	Only use non-sparking equipment.
P243	Take precautionary measures against static electricity.
P261	Avoid inhaling mist / vapor / aerosols.
P264	Wash hands thoroughly after handling.
P270	Do not eat, drink or smoke when using the product.
P271	Use only in open air or in a well-ventilated area.
P280	Use protective clothing/protective gloves/eye protection/face protection.
P312	In case of health problems, call the National Toxicological Information Center or a doctor.
P330	Rinse your mouth.
P363	Wash contaminated clothing before using it.
P301 + P312	Following ingestion: If you have any health problems, call the National Toxicology Information Center or a doctor immediately.
P302+P352	Following skin contact: Wash with plenty of soap and water.
P303+P361+P353	Following skin contact: Remove/take off all contaminated clothing. Immediately rinse the skin with water.
P304+P340	Following inhalation: Move the victim to fresh air and keep in a relaxed, motionless position that allows comfortable breathing.
P370+P378	In the event of fire: If it's safe, stop the leakage.
P403+P235	Keep in a well-ventilated place. Store in a cool place.
P501	Dispose of contents/container in accordance with local regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

Type of inspection	Inoculated culture	Result
Gram-positive	Staphylococcus aureus CCM 3953	blue-violet coca
Gram-negative	Escherichia coli CCM 3954	pink rods

Manufacturer
EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:
Obchodná ulica 7
946 51 Nesvady
Slovak Republic