

## Preparation Solution for the KATO Method (KATO)

In vitro diagnostic medical aid.

**Only for professional use.**

**Purpose:** The product is designed for impregnation of cellophane strips for KATO method for parasitological examination of stools.

**Package contents:** 1000 ml solution in a bottle

**Composition:** malachite green (C.I. 42000), glycerine, phenol,  
demineralized water

**Storage:** Store in a dark bottle below +25 °C.

**Usability:** Use within 30 days after opening the primary packaging.

**Procedure:** Use the solution according to internal procedures.

**Information for users:**

Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

**Warning:**



Warning

H226 Flammable liquid and vapor.  
H302 Harmful if swallowed  
H312 Harmful in contact with skin  
H332 Harmful if inhaled

P210 Keep away from heat / sparks / open flame / hot surfaces. Do not smoke.  
P233 Keep the container tightly closed.  
P240 Ground/fix the container and the filler.  
P241 Use an electrical / ventilation / lighting device designed for an explosive atmosphere.  
P242 Only use non-sparking equipment.  
P243 Take precautionary measures against static electricity.  
P261 Avoid inhaling gas / mist / vapor / aerosols.  
P264 Wash hands thoroughly after handling.  
P270 Do not eat, drink or smoke when using the product.  
P271 Use only in open air or in a well-ventilated area.  
P280 Use protective clothing/eye protection/protective gloves/face protection.  
P330 Rinse your mouth.  
P301 + P312 Following ingestion: If you have any health problems, call the National Toxicology Information Center or a doctor immediately.  
P302+P352 Following skin contact: Wash with plenty of soap and water.  
P363 Wash contaminated clothing before using it.  
P304+P340 Following inhalation: Move the victim to fresh air and keep in a relaxed, motionless position that allows comfortable breathing.  
P303+P361+P353 Following skin contact: Remove/take off all contaminated clothing. Immediately rinse the skin with water.  
P370+P378 In the event of fire: To extinguish use a foam/snow/water/powder fire extinguisher, depending on the materials in the immediate vicinity.  
P403+P235 Keep in a well-ventilated place. Store in a cool place.  
P501 Dispose of contents/container in accordance with local regulations.

**Instructions in case of a damaged packaging:**

In case of damage to the packaging, do not use and dispose of the product.

**Manufacturer**  
EnviroLab s.r.o.  
Galvaniho 7  
821 04 Bratislava  
Slovak Republic

**Place of business:**  
Obchodná ulica 7  
946 51 Nesvady  
Slovak Republic