

Potassium Nitrate 20% (KNO₃)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed for the preparation of disks to monitor the assimilation properties of the yeast (source of nitrogen).

Package contents: 1000 ml solution in a bottle

Composition: potassium nitrate
demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Warning

- H272 May contribute to the development of fire; oxidizing agent.
- P210 Keep away from heat / sparks / open flame. Do not smoke.
- P220 Keep away from flammable materials.
- P221 Adopt measures to prevent mixing with flammable materials.
- P280 Wear protective gloves, protective clothing, safety glasses.
- P370+P378 In the event of a fire, use a foam or snow fire extinguisher.
- P501 To extinguish use a snow, foam, water or powder fire extinguisher, depending on the materials stored in the immediate vicinity.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic