

## 0.5% Crystal Violet Solution (KRIVIO)

In vitro diagnostic medical aid.  
**Only for professional use.**

**Purpose:** The product is designed for Gram staining of microscopic preparations.

**Package contents:** 500 ml solution in a bottle

**Composition:** ethanol, crystalline violet (C.I. 42555), ammonium oxalate monohydrate  
demineralised water

**Storage:** Store in a dark bottle below +25 °C.

**Usability:** Use within 30 days after opening the primary packaging.

**Procedure:** Use the solution according to internal procedures.

**Information for users:** Do not use after the expiration date.  
Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

**Warning:**



Warning, Danger

H351	Suspected of causing cancer.
H302	Harmful if swallowed.
H312	Harmful in contact with skin.
H318	Causes serious eye damage.
H410	Very toxic to aquatic organisms, with long-lasting effects.
H225	Very flammable liquid and vapor.
P201	Read specific instructions before use.
P202	Do not use before you read and understand all safety precautions.
P264	Wash hands thoroughly after handling.
P270	Do not eat, drink or smoke when using the product.
P280	Use protective clothing/protective gloves/eye protection/face protection.
P273	Avoid release to the environment.
P308+P313	After exposure or suspected exposure: Seek medical help/care.
P301+P312	Following ingestion: If you have any health problems, call the National Toxicology Information Center or a doctor immediately.
P302+P352	Following skin contact: Wash with plenty of soap and water.
P363	Wash contaminated clothing before using it.
P305+P351+P338	Following eye contact: Gently rinse eyes with water for a few minutes. Remove contact lenses, if possible. Continue rinsing.
P391	Collect the released product.
P405	Keep locked up.
P501	Dispose of contents/container according to local regulations.

**Instructions in case of a damaged packaging:** In case of damage to the packaging, do not use and dispose of the product.

**Quality control:**

Type of inspection	Inoculated culture	Result
Gram-positive	Staphylococcus aureus CCM 3953	blue-violet coca
Gram-negative	Escherichia coli CCM 3954	pink rods

**Manufacturer**  
EnviroLab s.r.o.  
Galvaniho 7  
821 04 Bratislava  
Slovak Republic

**Place of business:**  
Obchodná ulica 7  
946 51 Nesvady  
Slovak Republic