
Lugol's Solution (LUGOL)

In vitro diagnostic medical aid.
Only for professional use.

Purpose: The product is designed for Gram staining of microscopic preparations.

Package contents: 1000 ml solution in a bottle

Composition: Potassium iodide
crystalline iodine
demineralized water

Storage: Store in a dark bottle below +25 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Quality control:

| Type of inspection | Inoculated culture | Result |
|--------------------|--------------------------------|-------------------|
| Gram-positive | Staphylococcus aureus CCM 3953 | blue-violet cocci |
| Gram-negative | Escherichia coli CCM 3954 | pink rods |

Manufacturer
EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:
Obchodná ulica 7
946 51 Nesvady
Slovak Republic