

Stock Solution for the MIF Method (MIF)

In vitro diagnostic medical aid.
Only for professional use.

Purpose: The product is designed to concentrate intestinal helminths and protozoa in stool specimens.

Package contents: 1000 ml solution in a bottle

Composition: Thimerosal, formaldehyde, glycerine, merthiolate
 demineralized water

Storage: Store in a dark bottle below +25 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.
 Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Warning

H317	May cause an allergic skin reaction
P261	Avoid inhaling smoke / gas / mist / vapor / aerosols.
P272	It is forbidden to take the contaminated work clothes away from the workplace.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352	Following skin contact: Wash with plenty of soap and water.
P333+P313	If skin irritation or rashes occur: Seek medical help/care.
P363	Wash contaminated clothing before using it.
P501	Dispose of contents/container in accordance with local regulations.

Instructions in case of a damaged packaging:

In case of damage to the packaging, do not use and dispose of the product.

Manufacturer
 EnviroLab s.r.o.
 Galvaniho 7
 821 04 Bratislava
 Slovak Republic

Place of business:
 Obchodná ulica 7
 946 51 Nesvady
 Slovak Republic