
1% Malachite Green Solution (MZ)

In vitro diagnostic medical aid.

Only for professional use.

Purpose:	The product is designed for staining of microscopic preparations.
Package contents:	1000 ml solution in a bottle
Composition:	malachite green (C.I. 42000), ethanol, demineralized water
Storage:	Store in a dark bottle below +25 °C.
Usability:	Use within 30 days after opening the primary packaging.
Procedure:	Use the solution according to internal procedures.
Information for users:	Do not use after the expiration date. Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Warning

H361 Suspected of harming the unborn child.

P201	Read specific regulations before use.
P202	Do not use before you read and understand all safety precautions
P281	Use the prescribed personal protective equipment
P308+P313	After exposure or suspected exposure: Seek medical help/care
P405	Keep locked up.
P501	Dispose of contents/container in accordance with local regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer
EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:
Obchodná ulica 7
946 51 Nesvady
Slovak Republic