
Phosphate Buffered Saline (PBS)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed for the dilution of samples in determination of phagocytic activity.

Package contents: 800 ml solution in a bottle

Composition: sodium chloride

disodium hydrogen phosphate

potassium dihydrogenphosphate

demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic