
Buffer Solution according to Kozak and Magrova (RKM1000)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed for examination of stools for helminths.

Package contents: 900 ml solution in a bottle

Composition: zinc sulphate
magnesium sulphate
glycerin
demineralized water

Storage: Store at a temperature of up to +25 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Danger

H302 Harmful if swallowed
H318 Causes serious eye damage.
H400 Very toxic to aquatic organisms.
H410 Very toxic to aquatic organisms with long-lasting effects.

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P264 Avoid inhaling dust/smoke/gas/mist/vapor/aerosols.
P270 Wear protective gloves / eye protection / face protection.
P273 Avoid release to the environment.
P305 + P351 + P338 FOLLOWING EYE CONTACT: Gently rinse eyes with water for a few minutes. Remove contact lenses, if possible. Continue rinsing.
P310 Call the NATIONAL TOXICOLOGY INFORMATION CENTER or a doctor immediately.
P301 + P314 FOLLOWING INGESTION: In the event of health problems, seek medical help/care.
P391 Collect the released product.
P501 Dispose of contents/container in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic