

Zinc Sulphate Solution 33% (SZ)

In vitro diagnostic medical aid.

Only for professional use.

Purpose: The product is designed to concentrate eggs of intestinal helminths and protozoa in stool specimens.

Package contents: 1000 ml solution in a bottle

Composition: zinc sulfate
demineralized water

Storage: Store in a cool place at a temperature of up to +25 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



H302	Harmful if swallowed.
H318	Causes serious eye damage.
H410	Very toxic to aquatic organisms, with long-lasting effects.
P273	Avoid release to the environment.
P290	Wear protective gloves / eye protection / face protection.
P305 + P351 + P338	FOLLOWING EYE CONTACT: Gently rinse eyes with water for a few minutes. Remove contact lenses, if possible. Continue rinsing.
P501	Dispose of the contents/container in an approved waste disposal facility.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic