

Tris-Borate Stock Solution (TBE)

In vitro diagnostic medical aid.
Only for professional use.

Purpose: The product is designed for the preparation of agarose gels and electrophoretic apparatus.

Package contents: 1000 ml solution in a bottle

Composition: tris(Hydroxymethyl)aminomethane

EDTA

boric acid

demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Warning:



Warning	H302	Harmful if swallowed
	H315	Irritating to skin
	H319	Causes serious eye irritation
	H335	May cause respiratory irritation

P261	Prevent inhalation of aerosols.
P264	Wash hands thoroughly after handling.
P270	Do not eat, drink or smoke when using the product.
P271	Use only in open air or in a well-ventilated area.
P280	Wear protective gloves/clothing/eye protection/face protection.
P301+P312	FOLLOWING INGESTION: In case of health problems, call the National Toxicological Information Center or a doctor.
P302+P352	FOLLOWING SKIN CONTACT: Wash with plenty of soap and water.
P304+P340	FOLLOWING INHALATION: Move the victim to fresh air and keep in a relaxed, motionless position that allows comfortable breathing.
P305+P351+P338	FOLLOWING EYE CONTACT: Gently rinse eyes with water for a few minutes. Remove contact lenses, if possible. Continue rinsing.
P312	In case of health problems, call the National Toxicological Information Center or a doctor immediately.
P330	Rinse your mouth.
P332+P313	If skin irritation occurs: Seek medical help/care.
P337+P313	If eye irritation persists: Seek medical help/care.
P362	Take off contaminated clothing and wash before using it.
P403+P233	Keep in a well-ventilated place. Keep the container tightly closed.
P501	Dispose of contents/container in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer
 EnviroLab s.r.o.
 Galvaniho 7
 821 04 Bratislava
 Slovak Republic

Place of business:
 Obchodná ulica 7
 946 51 Nesvady
 Slovak Republic